



SSTO6B211500

08/2024

Due to our policy of constant technological upgrading, technical specifications may be subject to change without prior notice.

According to the standards in force, in extra-EU areas the availability and specifications of some products and/or characteristics may vary. Please contact your local distributor for further information.



Making Your Life Better.

BU Medical Equipment
Sede Legale ed Amministrativa / Headquarters
CEFLA s.c.
Via Selice Provinciale, 23/a
40026 Imola (BO) - Italy
t. +39 0542 653111
f. +39 0542 653344
sternweber@sternweber.com
sternweber.com

Stabilimento / Plant
Via Bicocca, 14/C
40026 Imola (BO) - Italy
t. +39 0542 653441
f. +39 0542 653601



T7 T8 T9 STOOLS

COMFORT IS EVOLVING

STERNWEBER.COM

STOOLS
T7 T8 T9

COMFORT IS EVOLVING

2

THE PERFECT ADDITION TO ANY DENTAL PRACTICE

3





Seamless interaction with body and workplace

The new T7 stool interacts smoothly with the anatomy of the body, making the dentist's work easier by aiding proper posture at all times. High-quality materials and innovative design ensure treatment is always performed under optimal conditions. Its dimensions match the spaces available in a surgery, yet also leave full freedom of movement.

A comfortable seat - made of soft, easily sanitized material - helps dentists maintain proper positioning of legs and back during treatment.



Backrest
The backrest has been carefully sized to provide adequate lumbar support and minimise contact-generated heat.




Lumbar support
Independent adjustment of backrest height provides proper support for all dentists, whatever their build.



Tilt
The tilt angle can be adjusted to adapt the backrest to the dentist's own operating style and help decompress the spine when treatment is not being carried out.



Height adjustment
The T7 stool has an extensive height range, with independent adjustment of seat and back. A 10 cm extension is available to raise the seat further if required.




Perfect support for dentists

The T9 stool is the last word in ergonomics.

Thanks to revolutionary saddle-like geometry, the T9 combines a very small footprint with unmatched comfort. The vertical, open position of the legs makes it easier to approach the patient and maintain proper back posture. This reduces fatigue, even during prolonged treatment sessions.

Beneficial effects

Careful shaping ensures bodyweight is distributed optimally over the seat. This results in a more natural, relaxed working posture that keeps concentration and energy levels high.




Freedom of movement

A tilt mechanism ensures the 'saddle' supports the dentist's body movements. This dynamic seating solution favours a proper, active posture in every treatment scenario. The angle of movement is easily adjusted and can be locked in place.



Upholstery

Two T9 versions are available: with seamless upholstery for easier cleaning and with Memory Foam padding (optional), which offers exceptional comfort.



Optimal height

A special adjustment mechanism means the seat can easily be adapted to users of any height. The model includes optional accessories such as a 10 cm extension and/or footrest ring.

The stool that streamlines workflows

The T8 stool offers perfect support, even during long treatment sessions, while allowing full freedom of movement in every direction. Adjustment and personalisation mechanisms streamline adaptation to the surgery layout and enhance operator efficiency.



Rotating armrest and circular seat
The fully independent armrest can be rotated 360°. The seat is comfortable in any position.



Stability and practicality
The shaped armrest is horizontal and has no edges, providing better forearm support. Made of non-slip material that optimises grip, it is perfect for prolonged use.



Personalised support
During assembly the armrest can be positioned with a centred coupling, or off-centre to the right or left (to allow for prevalent use on one side or the other). Furthermore, the support can be adjusted at any time by rotating the armrest around the operator to the desired position.



	FOOTREST RING	EXTENSION 10 CM	BALANCE SYSTEM	STANDARD UPHOLSTERY	MEMORY FOAM PADDIANG
T7	○	○	-	●	-
T8	○	○	-	●	-
T9	○	○	●	●	○

● as standard * optional - unavailable




Footrest ring
An optional footrest ring providing comfortable, stable support is available on every model. A metallic finish ensures easy cleaning.



Extension
An optional 10 cm extension is available for all models, ensuring the seat height can be adapted to all user builds.



T7 STOOL



T8 STOOL



T9 STOOL



- | | | | | | |
|--|---|--|---|---|--|
| <p>STANDARD</p> <ul style="list-style-type: none"> 102 Atlantic blue 106 Mediterranean blue 113 Pacific blue 136 Indian blue 132 Blueberry violet 134 Japanese wisteria 135 Venetian red | <ul style="list-style-type: none"> 115 Scottish salmon 103 Nevada yellow 123 Polynesian green 101 Caribbean green 137 Bright silver 121 Anthracite grey 130 Graphite black | <ul style="list-style-type: none"> 142 Anatolian Hazelnut 141 Papyrus Beige 143 Arabian Gold 140 Brazilian Brown 144 Ruby Red | <p>MEMORY FOAM</p> <ul style="list-style-type: none"> 198 Atlantic blue 196 Mediterranean blue 183 Pacific blue 186 Indian blue 192 Blueberry violet 184 Japanese wisteria 194 Venetian red | <ul style="list-style-type: none"> 195 Scottish salmon 182 Nevada yellow 193 Polynesian green 197 Caribbean green 187 Bright silver 199 Anthracite grey 180 Graphite black | <ul style="list-style-type: none"> 172 Anatolian Hazelnut 171 Papyrus Beige 173 Arabian Gold 170 Brazilian Brown 174 Ruby Red |
|--|---|--|---|---|--|

Measurements in drawings are in millimetres